



Dual-Wavelength Retarder

Dual wavelength retarders can provide the same retardance at two wavelengths that are separated in wavelength by more than the span covered by an achromatic retarder. They can also provide different specified retardances at two different wavelengths.

Traditionally these retarders have been made using crystal quartz and are multiorder retarders at both wavelengths. Our dual wavelength retarders use polymers instead. They are usually much lower order and consequently have a slower change in retardance with angle of incidence as shown in the graph. On average the order is about 20% of that for a comparable quartz dual wavelength retarder.

Call for a quote on a custom coating on these normally uncoated retarders. The retardance tolerance is ± 0.01 waves at both wavelengths. Many custom combinations not listed in the catalog are available. Please call for a quote on your custom requirement. Standard unmounted sizes are 0.50 inches and 1.00 inches.



Key Features

• • •

Low order

Wide angular field

Broad wavelength coverage

Coated or Uncoated Available

Mounted or Unmounted Available

Waveplate Suite

• • •

Precision Retarder

Precision Achromatic Retarder

Precision Superachromatic Retarder

Dual-Wavelength Retarder

Wide Field Retarder

Liquid Crystal Variable Retarder

Polymer Film Retarder

Raptor Applied Polymer Retarder

Large Aperture Retarder

Bi-Crystalline Achromatic Retarder



SPECIFICATIONS

| | |
|---|--|
| Retarder Material | Birefringent Polymer |
| Substrate Material | N-BK7 |
| Retardance Accuracy | $\leq \lambda/100$ at both wavelengths |
| Transmitted Wavefront Distortion | $\leq \lambda/4$ |
| Beam Deviation | ≤ 1 arc-min |
| Reflectance (per surface) | $\sim 4\%$ at normal incidence |
| Storage Temperature | design dependent |
| Operating Temperature | design dependent |

Custom anti-reflection coatings to provide less than 0.5% reflection at both wavelengths are available. Please call your Meadowlark Optics sales engineer for a quote.

ORDERING INFORMATION

| Thickness in. (mm) | Dimensions in. (mm) | Part Number |
|---------------------------|-------------------------------|-------------------------------------|
| 0.14 (3.6 mm) | 0.50 (\varnothing 12.7 mm) | D R1 R2 – d – λ_1/λ_2 |
| 0.27 (6.9 mm) | 1.00 (\varnothing 25.4 mm) | D R1 R2 – d – λ_1/λ_2 |

R1,R2 = Q for quarter wave, H for half wave

d = Diameter (e.g. -100 is 1" OD)

λ_1, λ_2 = Design wavelengths. (e.g. -0488 is 488 nm)

DQH-200- λ_1, λ_2 => 2" Quarter Waveplate @ λ_1 , Half @ λ_2 .