

Imaging in Scattering or Turbid Media

Overview: Adaptive optics was first utilized to correct for aberrations that are introduced when imaging through atmospheric turbulence. In monochromatic imaging systems or laser communication systems wavefront correction is most easily accomplished by adding a liquid crystal spatial light modulator to the imaging system. By applying an equal and opposite phase to the SLM it is possible to restore diffraction limited images. In recent years, much of the research on atmospheric turbulence correction is translating to biology, where biological systems introduce scattering and turbidity. For example, SLMs can be used in STED microscopes for deep tissue imaging. In order to maintain the structure of the excitation and depletion sources, the aberrations that the sources will encounter when passing through the sample must be pre-corrected for. Similarly SLMs used in multi-photon imaging systems are used to pre-correct for scattering and aberrations the illumination will encounter when exciting deep tissue targets.




Figure 1 Hemphill, A. S., Tay, J. W., & Wang, L. V. (2016). Hybridized wavefront shaping for high-speed, high-efficiency focusing through dynamic diffusive media. *Journal of biomedical optics*, 21(12), 121502.

Critical requirements: For this market the SLM must offer high resolution, phase stability, and high speed switching. The SLM resolution determines the ability to correct for complex aberrations. High phase stability ensures temporally stable excitation which is important when imaging in scattering media with significant losses. High speed SLMs allow for real time adaptive optics.

Recommended References:

1. Tzang, O., Caravaca-Aguirre, A. M., Wagner, K., & Piestun, R. (2018, May). Adaptive Wave-front shaping in Linear and Nonlinear Complex Media. In *CLEO: QELS_Fundamental Science* (pp. FF3H-5). Optical Society of America.
2. Zhuang, B., Xu, C., Geng, Y., Zhao, G., Chen, H., He, Z., & Ren, L. (2018). An Early Study on Imaging 3D Objects Hidden Behind Highly Scattering Media: a Round-Trip Optical Transmission Matrix Method. *Applied Sciences*, 8(7), 1036.
3. Turtaev, S., Leite, I. T., Mitchell, K. J., Padgett, M. J., Phillips, D. B., & Čižmár, T. (2017). Comparison of nematic liquid-crystal and DMD based spatial light modulation in complex photonics. *Optics express*, 25(24), 29874-29884.
4. Knarr, S. H., Lum, D. J., Schneeloch, J., & Howell, J. C. Compressive direct imaging of a billion-dimensional optical phase space.

5. Duan, D. Y., & Xia, Y. J. (2018). Computational ghost imaging with nondegenerate wavelength light source. *arXiv preprint arXiv:1801.10045*.
6. Duan, D. Y., & Xia, Y. J. (2018). Computational ghost imaging with nonlocal quantum correlations. *arXiv preprint arXiv:1801.10045*.
7. Schneider, J., & Aegerter, C. M. (2018). Dynamic light sheet generation and fluorescence imaging behind turbid media. *Journal of the European Optical Society-Rapid Publications*, 14(1), 7.
8. Burgi, K. W., Marciniak, M. A., Nauyoks, S. E., & Oxley, M. E. (2017, August). Exploiting redundant phase information of a reflection matrix. In *Optical Trapping and Optical Micromanipulation XIV* (Vol. 10347, p. 103470K). International Society for Optics and Photonics.
9. Liu, Y., Ma, C., Shen, Y., Shi, J., & Wang, L. V. (2017). Focusing light inside dynamic scattering media with millisecond digital optical phase conjugation. *Optica*, 4(2), 280-288.
10. Hemphill, A. S., & Wang, L. V. (2016, March). Hybrid iterative wavefront shaping for high-speed focusing through scattering media. In *Adaptive Optics and Wavefront Control for Biological Systems II* (Vol. 9717, p. 97170V). International Society for Optics and Photonics.
11. Hemphill, A. S., Tay, J. W., & Wang, L. V. (2016). Hybridized wavefront shaping for high-speed, high-efficiency focusing through dynamic diffusive media. *Journal of biomedical optics*, 21(12), 121502.
12. Zhuang, B., Xu, C., Geng, Y., Zhao, G., Zhou, L., He, Z., & Ren, L. (2017, June). Imaging through a scattering medium based on spatial transmission matrix. In *European Conference on Biomedical Optics* (p. 104160D). Optical Society of America.
13. Zhuang, B., Xu, C., Geng, Y., Zhao, G., Chen, H., & Ren, L. (2018, November). Imaging through highly scattering media based on optical transmission matrix. In *Tenth International Conference on Information Optics and Photonics* (Vol. 10964, p. 109642F). International Society for Optics and Photonics.
14. Zhang, Y., Wu, C., Song, Y., Si, K., Zheng, Y., Hu, L., ... & Gong, W. (2019). Machine learning based adaptive optics for doughnut-shaped beam. *Optics Express*, 27(12), 16871-16881.
15. Hofer, M., & Brasselet, S. (2019). Manipulating the transmission matrix of scattering media for nonlinear imaging beyond the memory effect. *Optics letters*, 44(9), 2137-2140.
16. Burgi, K. W., Marciniak, M. A., Nauyoks, S. E., & Oxley, M. E. (2016, September). Matrix methods for reflective inverse diffusion. In *Reflection, Scattering, and Diffraction from Surfaces V* (Vol. 9961, p. 996100). International Society for Optics and Photonics.
17. Burgi, K., Marciniak, M., Oxley, M., & Nauyoks, S. (2017). Measuring the reflection matrix of a rough surface. *Applied Sciences*, 7(6), 568.
18. Premillieu, E., & Piestun, R. (2019). Measuring the Transmission Matrix of a scattering medium using Epi-Fluorescence light. *arXiv preprint arXiv:1909.07955*.
19. Thendiyammal, A., Osnabrugge, G., Knop, T., & Vellekoop, I. M. (2020). Model-based wavefront shaping microscopy. *arXiv preprint arXiv:2002.05279*.
20. Ják, P., Šiler, M., Ježek, J., Trägårdh, J., Zemánek, P., & Čížmár, T. (2019, February). Multimode fiber transmission matrix obtained with internal references. In *Adaptive Optics and Wavefront Control for Biological Systems V* (Vol. 10886, p. 1088610). International Society for Optics and Photonics.
21. Fleming, A., Conti, C., & Di Falco, A. (2018). Nonlinear transmission matrices of random optical media. *arXiv preprint arXiv:1809.07077*.
22. Boniface, A., Gusachenko, I., Dholakia, K., & Gigan, S. (2019). Rapid broadband characterization of scattering medium using hyperspectral imaging. *Optica*, 6(3), 274-279.
23. Ringuette, D., Sigal, I., Gad, R., & Levi, O. (2015). Reducing misfocus-related motion artefacts in laser speckle contrast imaging. *Biomedical optics express*, 6(1), 266-276.
24. Burgi, K., Ullom, J., Marciniak, M., & Oxley, M. (2016). Reflective inverse diffusion. *Applied Sciences*, 6(12), 370.
25. Zhuang, B., Xu, C., Geng, Y., Zhao, G., Chen, H., He, Z., ... & Ren, L. (2018). Round-trip imaging through scattering media based on optical transmission matrix. *Chinese Optics Letters*, 16(4), 041102.
26. Pozzi, P., Gandolfi, D., Porro, C. A., Bigiani, A., & Mapelli, J. (2020). Scattering Compensation for Deep Brain Microscopy: The Long Road to Get Proper Images. *Frontiers in Physics*, 8, 26.
27. Yin, X. L., Xia, Y. J., & Duan, D. Y. (2018). Theoretical and experimental study of the color of ghost imaging. *Optics express*, 26(15), 18944-18949.

28. Motz, A. M. A., Czerski, J., Adams, D. E., Durfee, C., Bartels, R., Field, J., ... & Squier, J. (2020). Two-dimensional random access multiphoton spatial frequency modulated imaging. *Optics Express*, 28(1), 405-424.
29. Zhao, M., Wang, H., & Tian, Z. (2019). Wavefront-shaping-based pattern regeneration through the scattering medium. *JOSA A*, 36(9), 1483-1487.
30. Escobet-Montalbán, A., Wijesinghe, P., Chen, M., & Dholakia, K. (2019, February). Wide-field multiphoton imaging with TRAFIX. In *Multiphoton Microscopy in the Biomedical Sciences XIX* (Vol. 10882, p. 108821G). International Society for Optics and Photonics.